

Cell Transfection Reagents

DESCRIPTION

The manipulation of gene expression is an important tool for researchers, yet primary cells are notoriously difficult to transfect. ScienCell Research Laboratories' suite of cationic polymer-based transfection reagents are designed and optimized for specific primary cell types. Obtain efficient transfection in just one day without the need to remove serum or antibiotics.

Products:

- [AstroFectagen](#)
- [EndoFectagen](#)
- [EpiFectagen](#)
- [EpiFectagen II](#)
- [FibroFectagen](#)
- [KeraFectagen](#)
- [MesenFectagen](#)
- [SMCFectagen](#)

