

BA-4/VT, Bucket adapters

DESCRIPTION

Bucket adapters for 4 vacutainer adapters of 2-9 ml (set of 4 pcs.). Included 8-9 ml [16 x 100] vacutainer adapters, BN-16/100-16.

SPECIFICATIONS

Tube holders material	POM Kocetal
Rotor type	Swing-out
Capacity	4
Tube's volume	2-9 ml
Max. speed	4200 rpm
Max. RCF	2960 × g



CAT. NUMBER

BS-010221-4K

BA-4/VT



LMC-56
Without rotor
Laboratory Centrifuge

LMC-56 is a modern benchtop, low-speed laboratory centrifuge designed for convenient sedimentation, centrifugation and collection of necessary samples. It provides operation with tubes, blood collection systems, gel cards, microtest and ...

[read more](#)



BR-4U
BS-010221-AK
Rotor

Bucket rotor without bucket adapters



BN-13/75-16
BS-010221-OK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/75-16 for vacutainers 2-5 ml [13x75]



BN-13/100-16
BS-010221-PK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-13/100-16 for vacutainers 4-8 ml [13x100]



BN-16/100-16
BS-010221-QK
adapter set

For bucket adapters BA-4/VT optional adapter set (16 pcs) - BN-16/100-16 for vacutainers 8-10 ml [16x100]



RR-U
BS-010208-UK
rack for rotors

Rack for rotors