

# PS R1, Precision Balance

## DESCRIPTION

### Communication interfaces

New PS.R1 balances, have pans in two possible dimensions: 128x128 mm or 195x195 mm balances with a smaller pan have a draft shield. They were equipped with the system of external mass adjustment.

PS.R1 balances feature several communication interfaces: 2 x RS 232, type A USB, type B USB and optional Wireless Connection.

### Databases

Information system of the PS R1 balances is based on 5 databases, which allows for several users to work with several products databases, and the registered weighing results can be subject to further analysis.

### ALIBI Memory

The R1 series balances feature ALIBI memory that is a warranty for safety and automatic recording of your measurements.

The ALIBI memory allows to record up to 100 000 weighing records.

[More info](#)

