

VP-8/5, block

DESCRIPTION

8 × 5 ml conical tubes

SPECIFICATIONS

Max. RPM

1200

Max. temperature

100°C



CAT. NUMBER

BS-010175-SK

VP-8/5



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)