

VP-4, block

DESCRIPTION

4 × 50 ml conical tubes

SPECIFICATIONS

Max. RPM	1000
Max. temperature	80°C
Block socket $\varnothing$	29.5 mm
Block socket depth	76.5 mm



CAT. NUMBER

BS-010175-GK

VP-4



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)