

VP-8/15, block

DESCRIPTION

8 × 15 ml conical tubes

SPECIFICATIONS

Max. RPM	1100
Max. temperature	80°C
Block socket $\varnothing$	17 mm
Block socket depth	80 mm



CAT. NUMBER

BS-010175-HK

VP-8/15



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)