

VP-32, block

DESCRIPTION

32 × 0.5 ml microtubes

SPECIFICATIONS

Max. RPM	1400
Max. temperature	100°C



CAT. NUMBER

BS-010175-JK	VP-32
--------------	-------



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)