

VP-CV-20, block

DESCRIPTION

20 × 10 mm cuvettes

SPECIFICATIONS

Max. RPM	1100
Max. temperature	80°C
Block socket $\varnothing$	□ 12.5 x 12.5 mm
Block socket depth	31 mm



CAT. NUMBER

BS-010175-IK

VP-CV-20



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)