

VP-CL-24, block

DESCRIPTION

24 × 3.6–4.5 ml cryotubes

SPECIFICATIONS

Max. RPM	1300
Max. temperature	100°C
Block socket $\varnothing$	12.5 mm
Block socket depth	30 mm



CAT. NUMBER

BS-010175-KK

VP-CL-24



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)