

VP-CS-24, block

DESCRIPTION

24 × 1–1.8 ml cryotubes

SPECIFICATIONS

Max. RPM	1400
Max. temperature	100°C
Block socket $\varnothing$	12.5 mm
Block socket depth	20 mm



CAT. NUMBER

BS-010175-LK

VP-CS-24



TS-100

Without thermoblock

Thermo-Shaker for microtubes and PCR plates

Thermo-Shaker **TS-100** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. Functions of heating (up to +100°C) and mixing can be performed both simultaneously and ...

[read more](#)