

VP-CV-20C, block

DESCRIPTION

20 × 10mm cuvettes

SPECIFICATIONS

Max. RPM	1100
Max. temperature	80°C
Block socket ø	□ 12.5 x 12.5 mm
Block socket depth	31 mm



CAT. NUMBER

BS-010176-IK

VP-CV-20C

ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)