

VP-CL-24C, block

DESCRIPTION

24 × 3.6–4.5 ml cryotubes

SPECIFICATIONS

Max. RPM	1300
Max. temperature	100°C
Block socket ø	12.5 mm
Block socket depth	30 mm



CAT. NUMBER

BS-010176-KK

VP-CL-24C

ACCESSORIES



TS-100C

Without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples down ...

[read more](#)



TS-100C Smart

Software included, without thermoblock

Thermo-Shaker with cooling for microtubes and PCR plates

Thermo-Shaker **TS-100C Smart** provides intensive mixing and temperature control of samples in microtest tubes or PCR plate. This model of Thermo-Shaker differs from TS-100 with a possibility of cooling samples ...

[read more](#)