

MSR-16

DESCRIPTION

Rotor for two 8 x 0.2 ml PCR strip rows or 16 individual 0.2 ml PCR tubes.

SPECIFICATIONS

Amount of 0.2 ml tubes	16
------------------------	----



CAT. NUMBER

BS-010218-AK	MSR-16
--------------	--------



Microspin 12 Plus

Including MSL-SC lid, adapters A-02, A-05
High-speed Mini-centrifuge

High-speed mini-centrifuge **Microspin 12 Plus** is a compact desktop centrifuge designed for biomedical laboratories.

Microspin 12 Plus is used for extracting RNA/DNA samples, sedimentation of biological components, biochemical and chemical

...

[read more](#)