

B1-L, block

DESCRIPTION

Mini water or beads bath, 55 x 100 x 73 mm

SPECIFICATIONS

Depth	73 mm
Capacity	max 0,4 L



CAT. NUMBER

BS-010418-WK

B1-L



CH 3-150
Without thermoblock
Heating and cooling thermostat

Combitherm-2 CH 3-150 is specially designed to thermostabilise materials at temperatures from -3°C to $+150^{\circ}\text{C}$ according to methods of analysis. To obtain useful functionality and decrease foot-print of instruments Combitherm-2 ...

[read more](#)