

SH-12x8, PCR strip tube holder set

DESCRIPTION

Magnetically fixed cover for easy fixing of various strip tube configuration.
PCR strip tube holders set of 2 pcs.

SPECIFICATIONS

Amount of tubes	96
-----------------	----



CAT. NUMBER

BS-010204-CK	SH-12x8, PCR strip tube holder set
--------------	------------------------------------



MSC-2P
BS-010204-A02
Minicentrifuge-Vortex for PCR plates

MSC-2P is a compact sized digital centrifuge intended to collect droplets, mix reagents and collect once more for improved PCR yield in subsequent analysis. The combination of spin-mix functions ensures ...

[read more](#)