

FS-1, Foot switch

DESCRIPTION

FS-1, Foot switch with USB-A connector. Wire length 2m.



CAT. NUMBER

BS-010177-AK

FS-1, Foot switch

ACCESSORIES



UIS-360

Without platform, protective screen, foot-switch
Universal Inoculation Spinner

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, ...

[read more](#)



FTA-U 4 L

Universal vacuum aspirator

FTA-U vacuum aspiration system is powerful, ergonomic, and durable benchtop vacuum system designed to improve liquid waste handling in demanding laboratory environments, including research, clinical diagnostics, pharmaceutical QC laboratories and ...

[read more](#)