


B-05G, block

DESCRIPTION

Block for one deep-well plate Greiner MASTERBLOCK® 96/0.5 ml - Cat.No. 786201, 786261



CAT. NUMBER

BS-010159-XK

B-05G block for one deep-well plate Greiner
MASTERBLOCK® 96/0.5 ml



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)