

Dust filter (2 pcs)

DESCRIPTION



CAT. NUMBER

BS-000000-S31

Dust filter (2 pcs)



UVC/T-AR
BS-040102-A13
PCR cabinet

The Biosan **UVC/T-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVC/T-M-AR
With built in socket and inlet for power cords
PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVT-B-AR
With built in socket and inlet for power cords
PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVT-S-AR
BS-040107-A07
PCR cabinet

The Biosan **UVT-S-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVR-M
BS-040105-AAA
UV Cleaner-Recirculator

- One 25W UV-lamp



UVR-Mi
BS-040110-AAA
UV Cleaner-Recirculator

- More powerful model with two 25W UV-lamps
- Timer feature
- LCD display