

UV lamp - 25W

DESCRIPTION



CAT. NUMBER

BS-000000-S09

UV lamp, 25W

ACCESSORIES



UVC/T-AR
 BS-040102-A13
 PCR cabinet

The Biosan **UVC/T-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVC/T-M-AR
 With built in socket and inlet for power cords
 PCR cabinet

The Biosan **UVC/T-M-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVT-B-AR
 With built in socket and inlet for power cords
 PCR cabinet

The Biosan **UVT-B-AR** is a benchtop DNA/RNA UV-cleaner workstation designed to create a clean, contamination-free environment for PCR setup and other nucleic acid-related workflows.

This model features a dual UV ...

[read more](#)



UVR-M
 BS-040105-AAA
 UV Cleaner-Recirculator

- One 25W UV-lamp



UVR-Mi
 BS-040110-AAA
 UV Cleaner-Recirculator

- More powerful model with two 25W UV-lamps
- Timer feature
- LCD display