

HOSE10, 10mm hose kit for the ORCA

DESCRIPTION

Includes 2x 2m hose (ID 10mm, 14.5mm OD), 2x 2m insulation, 4x hose clips and 2x hose barb (2x 1/2" male to 10mm barb (1/2" to 3/8"))



CAT. NUMBER

HOSE10

10mm hose kit for the ORCA



ORCA

ORCA EU

Closed loop refrigerated circulator

A new sustainability focused, energy-efficient recirculating chiller designed to replace tap-water-based cooling in laboratories.

Key Features

- Perfectly designed to replace tap water-based cooling solutions in labs and slash utility spend ...

[read more](#)