

## HOSE13, 13mm hose kit for the ORCA

### DESCRIPTION

Includes 2x 2m hose (ID 13mm, 19mm OD), 2x 2m insulation, 4x hose clips and 2x hose barbs (2x 1/2" male to 13mm barb (1/2"))



### CAT. NUMBER

HOSE13

13mm hose kit for the ORCA



**ORCA**

ORCA EU

Closed loop refrigerated circulator

A new sustainability focused, energy-efficient recirculating chiller designed to replace tap-water-based cooling in laboratories.

Key Features

- Perfectly designed to replace tap water-based cooling solutions in labs and slash utility spend ...

[read more](#)