

P-90, round platform for 90 / 105 mm Petri dishes

DESCRIPTION

The P-90 interchangeable round platform for the UIS-360 is designed for the convenient placement of standard 90 mm or 105 mm Petri dishes and for performing inoculations. The silicone mat and adapter ring are not included in the kit; the mat and adapter ring from the standard universal platform are intended for use.



CAT. NUMBER

BS-010177-CK

P-90



UIS-360
BS-010177-A01
Universal Inoculation Spinner

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, ...

[read more](#)