

P-90, round platform for 90 / 105 mm Petri dishes

DESCRIPTION

The P-90 interchangeable round platform for the UIS-360 is designed for the convenient placement of standard 90 mm or 105 mm Petri dishes and for performing inoculations.



CAT. NUMBER

BS-010177-CK

P-90



UIS-360

Without platform, protective screen, foot-switch
Universal Inoculation Spinner

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, ...

[read more](#)