

TSP-24/1.5, block

DESCRIPTION

24 x 1.5 ml microtubes

SPECIFICATIONS

Amount of 1.5 ml tubes	24
Max. RPM	2500
Max. temperature	100°C



CAT. NUMBER

BS-010178-KK

TSP-24/1.5



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)