

## TSP-24/2.0, block

### DESCRIPTION

24 × 2 ml microtubes

### CAT. NUMBER

BS-010178-JK

TSP-24/2.0

### SPECIFICATIONS

|                      |       |
|----------------------|-------|
| Amount of 2 ml tubes | 24    |
| Max. RPM             | 2500  |
| Max. temperature     | 100°C |

### ACCESSORIES



### TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)