

TSP-96/PCR, block

DESCRIPTION

96-well unskirted or semi-skirted microplate (0.2 ml) for PCR or 12 × 8 - 0.2ml strips or 96 tubes of 0.2 ml.

CAT. NUMBER

BS-010178-HK

TSP-96/PCR

SPECIFICATIONS

96-well unskirted or semi-skirted microplate (0.2 ml)	1
Max. RPM	2500
Max. temperature	100°C

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)