

TSP-20/0.5+12/1.5, block

DESCRIPTION

20 × 0.5 ml + 12 × 1.5 ml microtubes

SPECIFICATIONS

Amount of 0.5 ml tubes	20
Amount of 1.5 ml tubes	12
Max. RPM	2500
Max. temperature	100°C



CAT. NUMBER

BS-010178-MK

TSP-20/0.5+12/1.5



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)