

TSP-20/0.2+12/1.5, block

DESCRIPTION

20 × 0.2 ml + 12 × 1.5 ml microtubes

CAT. NUMBER

BS-010178-LK

TSP-20/0.2+12/1.5

SPECIFICATIONS

Amount of 0.2 ml tubes	20
Amount of 1.5 ml tubes	12
Max. RPM	2500
Max. temperature	100°C

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)