

TSP-32/0.5, block

DESCRIPTION

32 × 0.5 ml microtubes

CAT. NUMBER

BS-010178-EK

TSP-32/0.5

SPECIFICATIONS

Amount of 0.5 ml tubes	32
Max. RPM	2500
Max. temperature	100°C

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)