

TSP-8/5, block

DESCRIPTION

8 × 5 ml conical tubes

SPECIFICATIONS

Max. RPM

1200

Max. temperature

100°C



CAT. NUMBER

BS-010178-IK

TSP-8/5



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)