

TSP-8/15, block

DESCRIPTION

8 × 15 ml conical tubes

CAT. NUMBER

BS-010178-FK

TSP-8/15

SPECIFICATIONS

Max. RPM	1100
Max. temperature	80°C
Block socket $\varnothing$	17 mm
Block socket depth	80 mm

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)