

TSP-4/50, block

DESCRIPTION

4 × 50 ml conical tubes

SPECIFICATIONS

Max. RPM	1000
Max. temperature	80°C
Block socket $\varnothing$	29.5 mm
Block socket depth	76.5 mm



CAT. NUMBER

BS-010178-GK

TSP-4/50



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)