


TSP-20/CV10, block

DESCRIPTION

20 × 10 mm cuvettes

SPECIFICATIONS

Max. RPM	1100
Max. temperature	80°C
Block socket $\varnothing$	□ 12.5 x 12.5 mm
Block socket depth	31 mm



CAT. NUMBER

BS-010178-DK

TSP-20/CV10



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)