

TSP-20/CV10, block

DESCRIPTION

20 × 10mm cuvettes

CAT. NUMBER

BS-010178-DK

TSP-20/CV10

SPECIFICATIONS

Max. RPM	1100
----------	------

Max. temperature	80°C
------------------	------

Block socket $\varnothing$	□ 12.5 x 12.5 mm
----------------------------	------------------

Block socket depth	31 mm
--------------------	-------

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)