

TSP-24/CRY04, block

DESCRIPTION

24 × 3.6–4.5 ml cryotubes

CAT. NUMBER

BS-010178-CK

TSP-24/CRY04

SPECIFICATIONS

Max. RPM	1300
Max. temperature	100°C
Block socket $\varnothing$	12.5 mm
Block socket depth	30 mm

ACCESSORIES



TS-100C Pro

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)