

TSP-24/CRYO1.8, block

DESCRIPTION

24 × 1–1.8 ml cryotubes

SPECIFICATIONS

Max. RPM	1400
Max. temperature	100°C
Block socket $\varnothing$	12.5 mm
Block socket depth	20 mm



CAT. NUMBER

BS-010178-BK

TSP-24/CRYO1.8



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)