

Heated lid for TS-100C Pro

DESCRIPTION

Detachable independent heated cover, for thermal and condensation control.

Compatible with blocks:

1. [TSP-20/0.5+12/1.5](#)
2. [TSP-20/0.2+12/1.5](#)
3. [TSP-32/0.5](#)
4. [TSP-24/1.5](#)
5. [TSP-24/2.0](#)
6. [TSP-20/CV10](#)
7. [TSP-96/PCR](#)



CAT. NUMBER

BS-010178-NK

Heated lid



TS-100C Pro

BS-010178-A01

Programmable Thermo-Shaker with Heating and Cooling

The **TS-100C Pro** is a programmable laboratory thermo-shaker that combines precise heating and cooling, intensive mixing and interchangeable thermoblocks in one compact instrument.

Designed for repeatable sample-preparation workflows, it allows ...

[read more](#)