

U-150, universal platform

DESCRIPTION

The U-150 interchangeable universal platform fits both square (125 × 125 mm) and round Petri dishes (Ø90, Ø105, Ø150).



CAT. NUMBER

BS-010177-BK

U-150



UIS-360

Without platform, protective screen, foot-switch
Universal Inoculation Spinner

The **UIS-360** Universal Inoculation Spinner is a versatile laboratory device designed to streamline microbial plating and ensure consistent results. It provides smooth, uniform rotation for spreading samples across agar surfaces, ...

[read more](#)