


B-2MB, block

DESCRIPTION

Block for one deep-well plate KingFisher™ Apex Deepwell Plate - Cat.No. 95040450B



CAT. NUMBER

BS-010159-4K

Block for one deep-well plate KingFisher™ Apex
Deepwell Plate



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)