


B-2FN, block

DESCRIPTION

Block for one deep-well plate NEST Biotechnology® 96/2200 µl, I-Shaped - Cat.No 503761



CAT. NUMBER

BS-010159-3K

Block for one deep-well plate NEST
Biotechnology® 96/2200 µl, I-Shaped



TS-DW

Without thermoblock

Thermo-Shaker for deep well plates

TS-DW Thermo-Shaker is designed for shaking and incubating deep well plates.

A multisystem principle, used in the design of the Thermo-Shaker, allows operating it as 3 independent devices: Incubator, Plate

...

[read more](#)